

# DPC 2026 - Freiburg

Distangling contaminants from true intratumoural microbial signals in pancreatic ductal adenocarcinoma – A benchmarking approach


Linh Dang

2026-03-07

# Contamination in Microbiome Study

Cell

## Tumor Microbiome Diversity and Composition Influence Pancreatic Cancer Outcomes




Riquelme et al. Cell 2019

Key taxa: *Streptomyces*, *Pseudoxanthomonas*, *Saccharopolyspora*

Article | [Open access](#) | Published: 31 August 2021

## Tumor microbiome contributes to an aggressive phenotype in the basal-like subtype of pancreatic cancer

[Wei Guo](#), [Yuchao Zhang](#), [Shiwei Guo](#), [Zi Mei](#), [Huiping Liao](#), [Hang Dong](#), [Kai Wu](#), [Haocheng Ye](#), [Yuhang Zhang](#), [Yufei Zhu](#), [Jingyu Lang](#), [Lindian Hu](#) , [Gang Jin](#)  & [Xiangyin Kong](#) 




Guo Communication biology 2021

Key taxa: *Acinetobacter*, *Pseudomonas*, *Sphingopyxis*


## Section 1

### Negative Control Survey in the Lab


# Negative Sample Collection




- In total 203 negative control samples over 4 years




# Negative Sample Collection




- In total 203 negative control samples over 4 years
  - 113 buffer controls,




# Negative Sample Collection




- In total 203 negative control samples over 4 years
  - 113 buffer controls,
  - 84 paraffin controls,



# Negative Sample Collection




- In total 203 negative control samples over 4 years
  - 113 buffer controls,
  - 84 paraffin controls,
  - 6 PCR control.



# Bacterial Composition


## Composition of NCT Samples

Order by sum







# Heatmap from Negative Controls



## Section 2


# Intratumor Microbiome in PDAC and Assessment of

# Decontamination Strategies




Davis et al. 2018 Microbiome

## SCRuB




Austin et al. 2023 Nature Biotechnology



Nejman et al. 2020 Science


# Assessment of Decontamination Approach

- **Yield:**  
$$\frac{\# \text{ putative true taxa } \notin \text{NCT}}{\# \text{observed species}}$$




# Assessment of Decontamination Approach

- **Yield:**  
$$\frac{\# \text{ putative true taxa} \notin NCT}{\# \text{ observed species}}$$
- **Purity:**  
$$\frac{\# \text{ putative true taxa} \notin NCT}{\# \text{ putative contaminants}}$$




# Assessment of Decontamination Approach

- **Yield:**  
$$\frac{\# \text{ putative true taxa} \notin NCT}{\# \text{ observed species}}$$
- **Purity:**  
$$\frac{\# \text{ putative true taxa} \notin NCT}{\# \text{ putative contaminants}}$$
- **Composite Score =**  
$$\text{Yield} \times \text{Purity}$$




# Integrate Replica

## Utilize Replica for Assessment




- 10 PDAC sample, each with 2 replica

Paired Wilcoxon test  
The lower the better




# Integrate Replica

## Utilize Replica for Assessment




- 10 PDAC sample, each with 2 replica
- Aichson distance with paired Wilcox test.

Paired Wilcox test  
The lower the better








# Overlap - Microbiome profile of PDAC in mice




## Overlap - Microbiome profile of PDAC in mice




## Overlap - Microbiome profile of PDAC in mice



# Core microbiome after Nj approach



# Core microbiome after Nj approach



## Decontamination Methods

- Negative control samples are essential to avoid spurious taxa.

## Intratumor microbiome in PDAC

## Decontamination Methods

- Negative control samples are essential to avoid spurious taxa.
- Each batch has its own contaminant profile.

## Intratumor microbiome in PDAC

## Decontamination Methods

- Negative control samples are essential to avoid spurious taxa.
- Each batch has its own contaminant profile.
- In our study, thoroughly remove contaminant such as Nejman et al. procedure outperforms others.

## Intratumor microbiome in PDAC



## Decontamination Methods

- Negative control samples are essential to avoid spurious taxa.
- Each batch has its own contaminant profile.
- In our study, thoroughly remove contaminant such as Nejman et al. procedure outperforms others.

## Intratumor microbiome in PDAC

- Re-sequencing and various decontamination methods with proper negative controls could yield a reliable set of intratumor microbiota.

## Decontamination Methods

- Negative control samples are essential to avoid spurious taxa.
- Each batch has its own contaminant profile.
- In our study, thoroughly remove contaminant such as Nejman et al. procedure outperforms others.

## Intratumor microbiome in PDAC

- Re-sequencing and various decontamination methods with proper negative controls could yield a reliable set of intratumor microbiota.
- Almost identified taxa are well-known gut microbiota.

# Acknowledgements

UNIVERSITÄTSMEDIZIN  
GÖTTINGEN **UMG**

## AG Neesse

Christoph Ammer-Herrmenau

Jacob Hamm

Jutta Blumberg

Ulrike Wegener

Sören Buchholz

Catharina Lingsens

Cathleen Heuer

Areej Shahbaz

Richard Meier

Rahel Scholz

Louisa Eskelson

Helena Lescow

Charlotte Vogel

Linh Dang

Isabel Toussaint

Johannes Richter

Annemarie Dresbach

Leoni Florschütz



## UMG

AG Beissbarth (Bioinformatik)

AG Ströbel (Pathologie)

AG Schneider (Chirurgie)

AG Papantonis (Pathologie)

AG Sax (Medizinische Informatik)

## Externe Kooperationspartner

Johannes Zimmermann, Jena

Gabriel Nussbaum, Jerusalem

Bernd Heinrich, Hannover

Magdalena Huber, Marburg

Matthias Gaida, Mainz

Paco Real, CNIO Madrid

Nuria Malats, CNIO Madrid

Evgenij Shumilov, Münster

Ulrike Bacher, Bern



**P-MAPS Consortium**

## Section 3

### Backup Slides

# Core microbiome - All decontamination method - phylum: family

